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


ENGLISH

INTRODUCTION


- Super large LCD screen display, it is clear and easy to be read
- Fast weighing, the speed of weighing is from 10 times to 30 or 40 times faster than general mechanical balance
- It is easy to operate, weigh directly
- High intelligence: within the weighing range you can tare, zero clearing and accumulated weighing; over loading display and fault alarm
- Built-in hanger for below-balance weighing (Optional choice)
- In order to print the weighing data off through the printer or input the data to computer or other external devices, the comparator has RS-232 interface.

THE STRUCTURE OF BALANCE

The external view of balance:

Keyboard and display:

No	Display Sign	Meaning
1	g	gram
2	ct	carat
3	lb	pound
4	oz	ounce
5	•••••	piece counting mode
6	%	percentage mode
7	O	stable indicator
8	--End--	calibration is ended
9	Err--0	when calibration is "0" error
10	Err--1	
11	E	over the capacity
12	-E	under the readability
13	balance is processing data

CARE AND MAINTENANCE, ENVIRONMENTAL REQUIREMENTS

- Put the balance on a stable and fixed operating platform
- The operating platform should not be affected by excessive air flow caused by open window or a door
- The operating platform should not be affected by vibration. Four corners of the room are affected by vibration very little, perfect for setting up the operating platform; do not put the balance in direct sunshine
- Do not put the balance near heater or air conditioner
- Do not put the balance in areas where the danger of explosion exists
- Do not put the balance under raining or use water to clean, if there is water on the balance you should use dry cloth to dry, when the balance are not working properly, you should send it to repair
- Do not beat or strike or press the balance very hard
- Do not let tiny little object go inside the balance
- The best operating temperature $20^{\circ}\text{C} \pm 5^{\circ}\text{C}$
- Do not use a sharp or rough object (such as pencil or ball pen etc.) to poke the keys
- Do not let any object fall on the pan, otherwise the weighing system will be damaged
- Do not expose the balance in high temperature or mill dust environment for a long time
- Keep the balance clean and dry

Cleaning

- Unplug the AC adaptor before cleaning
- Do not use any aggressive cleaning agent (such as solvent). Use a piece of wet smooth and soft cloth with some mild detergent (such as soap) to clean
- Make sure no liquid or other tiny insects enter into the balance
-

OPERATION, CALIBRATION

Before using the balance, it should be placed on a stable and fixed operating platform. Level the balance by turning the adjusting feet, checking the level indicator on the balance, until the bubble appears in the center of the circle.

CALIBRATION:

- Plug in and pre-warm the balance at least 10 minutes
- Press CAL key to calibrate
- "C5000" is displayed
- Put the external calibration weight (5kg) on the pan carefully
- Press TARE key (**↕T**), wait the balance to collect data, after collecting data the value of external calibration weight is displayed (C5000 is flashing), "-End" is displayed, 5000.0g should be displayed; if the displayed value is not correct, you need to calibrate the balance again; take the weight off and repeat the procedure until displayed value is correct
- Take the weight off

WEIGHING UNIT EXCHANGE

There are 4 types of weighing units and 2 types of weighing modes available to choose

Press  key to change from one to another.
Order:

- Unit : Gram
- Unit : Carat
- Unit : Ounce
- Unit : Pound
- Mode : piece counting
- Mode : percentage

İşaret	İsim	Çevrim
g	Gram	1g
oz	Ounce	28.3495231g
dwt	Penny weight	1.55517384g
GN	Grain	0.06479891g
lb	Pound	453.7205g
ct	Carat	0.2g
kg	kilogram	1000g
tr	Tola (India)	11.6638038g
tb	Tael (Taiwan)	37.5g
tl	Tael (HK)	37.794g

WEIGHING

- No-load on the pan
- Press TARE key ()
- "0.0" is displayed
- Put the object on the pan
- When stable indicator "O" is displayed, read weighing value
- Remove the object, "0.0g" is displayed

WEIGHING WITH CONTAINER

- Put the container on the pan first
- The weight of container is displayed
- Press TARE key () to take the weight of container off, "0.0g" is displayed
- Put the object in the container
- when stable indicator "O" is displayed, read the value



PIECE COUNTING MODE

This mode is for counting standard objects (the objects have same weight):

- Press TARE key , "0.0g" is displayed
- Press MODE key  to chose piece counting mode ()
- Put 10 samples on the pan
- Press CAL key, "SEL..." is displayed first, then "10" (10 pieces) is displayed
- Remove all the samples, "0" (0 piece) is displayed
- Put some objects on the pan
- The quantity of objects is displayed

You can change number of samples among 10, 20, 50, 100 and 1000 by changing parameter setting C1-0, C1-1, C1-2, C1-3, C1-4

PERCENTAGE WEIGHING MODE %

- Press TARE key , "0.0g" is displayed
- Press MODE key  to chose percentage weighing mode (%)
- Put the reference weight or sample on the pan
- Press CAL key, "100.00" is displayed
- Remove reference weight or sample, "0.00" (0.00%) is displayed
- Put the contrast weight or sample on the pan
- The percentage of the contrast weight and reference weight is displayed

PARAMETER SETTING

PARAMETER SETTING TABLE:

Function	Setting	Meaning
C1 : Set the basic number of samples for piece counting	C1-0	10
	C1-1	20
	C1-2	50
	C1-3	100
	C1-4	1000
C2 : Background light control	C2-0	Auto
	C2-1	Always On
	C2-2	Always Off
C3 : "0" point minimum display value	C3-0	No "0" point tracking status
	C3-1	1d
	C3-2	2d
	C3-3	3d
	C3-4	4d
	C3-5	5d
C4 : Serial baud rate	C4-0	2400bps
	C4-1	1200bps
	C4-2	4800bps
	C4-3	9600bps
C5 : Data output mode	C5-0	On zero stable
	C5-1	On stable
	C5-2	On command or print key
	C5-3	Continuous
C6 : Key tone	C6-0	Yes
	C6-1	No
C7: Anti-interference	C7-0	Level 0
	C7-1	Level 1
	C7-2	Level 2

Factory default settings are C3-3, C5-2, C7-1, the other settings are 0

Parameter Setting Adjusting Method:

- Press  key more than 3 second
All LCD segments are displayed.
- Release and press again  key.
Display "C1-0" or similar.
- Change left number with  key, right number with  key.
- After setting all the parameters, press TARE key  to save the settings and go back to weighing status.

RS232 INTERFACE

- Baud Rate : (See table C4)
- Start bit : 1
- Data bits : 8
- Stop bit : 1
- Parity : No

DATA OUTPUT MODE:

See Table C5 parameter setting

- On zero stable output: "0.0000g" is displayed, put object on the balance, When displayed value is stable, a group of data will be transferred out
- Stable output: when weighing value is stable (stable indicator displayed), a group of data will be transferred out
- Command control output: PRINT key  and command P<CR> <LF> are print control commands, when the balance receive the command, a group of data will be transferred out
- Continuous output: every 0.3 sec a group of data will be transferred out
- PRINT key output mode: only when you press PRINT key , a group of data will be transferred out. Parameter setting is C5-2

Data output format:

FXXXXXX. XKKK<CR> <LF>

F : sign bit "+" or "-"

X : 0—9 weighing data

. : decimal point

K : Reserving three sign of weighing unit, right align, if there is not three sign space will fill in the room

<CR> : ENTER

<LF> : LINEFEED

For example, if the value is +5000.0g, after the transmission it will be +0005000.0g<CR><LF>

When it is unstable output, KKK are space; when it is stable output, there will be information sent out

Command control

This type of control mode can be used to control balance's [On/Off] key, TARE key, weighing unit and weighing mode's transformation, from external devices' commands.

- ① <O> <CR> <LF> [On/Off] command. It has the same function as comparator's [On/Off] key (4F 0D 0A);
- ② <T> <CR> <LF> TARE command. It has the same function as comparator's TARE key (54 0D 0A);
- ③ <C> <CR> <LF> CAL command. It has the same function as comparator's CAL key (43 0D 0A);
- ④ <M> <CR> <LF> Mode transformation command. It has the same function as comparator's MODE key (4D 0D 0A);
- ⑤ <P> <CR> <LF> Print command. It has the same function as balance's PRINT key. As long as the balance receive this command, a group of data will be transferred out (50 0D 0A).

Note: every time the balance receive a right command, it will respond the command and send the command back to external devices to indicate that the balance has received a right command from external devices.

Connection of balance and external device

Balance	Computer	Balance	Series interface printer
9 pins	9 holes	9 pins	25 pins
2	——— 2	2	——— 2
3	——— 3	3	——— 3
5	——— 5	5	——— 7

NOTE: the length of data line should not be more than 15 meters

TECHNICAL SPECIFICATIONS

Model	THT6K11	THT6102
Capacity	6100g	610g
Calibration weight	5000g	500g
Readability	0.1g	0.01g
Repeatability	±0.2g	±0.02g
Linearity	±0.2g	±0.02g
Stabilization time	about 2s	about 3s
Dimension	275x230x85 (L x W x H)mm	
Pan size	φ125	160X160X8mm
Weight	1.6kg	
Power	220V	50/60Hz

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SANAYİ VE TİCARET LTD.ŞTİ
MÜH. SAN. TIC. LTD. ŞTİ.**